



## Single-Use Cultivation in a Classical Format:

Reflecting the demand for single-use manufacturing by implementing stirred single-use bioreactors based on classical standards

Dr. Andreas Kocourek, Sartorius Stedim Biotech, Application Specialist, Fermentation Australia ([andreas.kocourek@sartorius-stedim.com](mailto:andreas.kocourek@sartorius-stedim.com))

## Overview

- Benefits & challenges of single-use bioreactors
- Single-use bioreactor with "classical" design, such as BIOSTAT CultiBag STR
- Some technical details
  1. Vessel design
  2. Bag design
  3. Impeller design
  4. Aeration/gassing design
- Case studies (mammalian cell culture)
- Summary



## Single-use bioreactors: Benefits



### Reduction:

- Contamination risk
- Time for set-up
- Space requirement
- Plant validation
- Investment
- Energy, WFI,...

### Increase:

- Flexibility

### Remove the need for:

- Cleaning & Sterilization

⇒ ↑ Process efficiency & ↓ Process Cost

# Impact of Single-Use to a manufacturing process: Case Study

## Stainless Steel Scenario

<b>Throughput</b>	Batches per year	<b>204 #/yr</b>		
	Annual product mass	<b>419.6 kg/yr</b>		
<b>Capital</b>	Total	<b>115,685,482 \$</b>		
	Per Litre of bioreactor	<b>9,640 \$/L</b>		
<b>CoG</b>	Cost per gram	<b>152.2 \$/g</b>		
<b>L/Batch</b>	<b>PW</b>	<b>WFI</b>	<b>Total</b>	<b>%</b>
Process	0	11,474	11,474	11%
Cleaning	55,869	37,191	93,060	89%
<b>Total</b>	<b>55,869</b>	<b>48,665</b>	<b>104,534</b>	<b>100%</b>

## Disposable Scenario

<b>Throughput</b>	Batches per year	<b>226 #/yr</b>		
	Annual product mass	<b>464.9 kg/yr</b>		
<b>Capital</b>	Total	<b>62,286,851 \$</b>		
	Per Litre of bioreactor	<b>5,191 \$/L</b>		
<b>CoG</b>	Cost per gram	<b>119.1 \$/g</b>		
<b>L/Batch</b>	<b>PW</b>	<b>WFI</b>	<b>Total</b>	<b>%</b>
Process	0	11,474	11,474	85%
Cleaning	1,600	450	2,050	15%
<b>Total</b>	<b>1,600</b>	<b>11,924</b>	<b>13,524</b>	<b>100%</b>



<b>Scenario</b>	<b>CoG \$/g</b>	<b>Capital \$million</b>
Stainless Steel	152.2	115.7
Disposable	119.1	62.3
<b>Savings</b>	<b>21.7%</b>	<b>46.2%</b>

Source: Biopharm Services

## Single-use bioreactors: Challenges for upscaling



### Challenges during upscaling

- Very complex process
- some parameters proved their usefulness are:
  - same geometry/ratios
  - impeller design
  - tip speed
  - Power input per volume (P/V)
  - Gassing strategy
  - kLa
- Knowledge and experience are nearly as important as the scaling-up parameters

### BUT


- Variations between reusable and currently existing single-use systems
- Single-use not existing for all volumes

# Single-use bioreactors: Difficulties to overcome

## Challenges in transfer due to:

- Change of Mixing Technology & behaviour
  - Wave induced mixing
  - Non-classical stirrer position
  - Non-classical impellers
- Gassing strategies
  - Only surface aeration
  - Non-classical spargers

⇒ Need for improvements in comparability reusable and disposable processes



From R. Eibl and D. Eibl; Adv Biochem Engin/Biotechnol (2009) 112: 183–207

# Single-use bioreactors with a "classical" design

## Case Study: BIOSTAT® CultiBag STR



## Requirements for „classical“ single-use bioreactors: Scaling up towards the new modern production chain



12,5-50 L

50-200 L

250-1000 L

# Requirements for „classical“ single-use bioreactors: Scaling up towards the new modern production chain



12,5-50 L



50-200 L



250-1000 L

## Requirements for „classical“ single-use bioreactors: control depending on application



### Flexibility by Application:

- For R&D, Seed, Pilot and Production
  - Basic & advanced version
  - Heating only or Heating/Cooling
  - Single & Twin set-up
  - Large variety of volumes:
    - » First systems 50L & 200L
    - » scalable from R&D to large scale production

## Requirements for „classical“ single-use bioreactors: advanced control with safe data-storage



# 1. Vessel design

## Classical requirements



### Single-Use Bioreactor Design: Bag Holder

- Rounded bottom as classical vessels
- Opening for harvesting port at lowest point
- Side openings for installation of probes
- Filter holders with filter heater for increase process safety



# 1. Vessel design

## Easy handling and time-saving




### BIOSTAT® STR Holder Design:

- Disconnectable from control tower for transfer to DSP and connection of another bag holder skid
  - ⇒ no time loss due to bag preparation, harvesting, further process steps
  - ⇒ very low operational downtime of equipment
- Double door for easy installation bag

# 1. Vessel Design Temperature Control

Heating blankets (heating only)



Double wall (heating and cooling)




# 1. Vessel Design


## Temperature Control with heating blanket (only heating)

**History Plot**  
 Aufheizen 4,5°C auf 37°C  
 CultiBag STR Temptest 2-2(Running) Selection :31.03.2009 14:50:14 - 01.04.2009 00:50:14

TEMP\_1.Ctrl\_output;Db 1.0 %  
 TEMP\_1.Setpoint;Db 0.50 °C  
 TEMP\_1.Value;Db 0.50 °C




# 1. New development: Vessel Design Temperature Control with double wall (heating & cooling)



## 2. Bag Design



## 2. Improved Bag Design: sterility, safety & longterm stability Flexibility in plastic



### CultiBag STR Design:

- Multilayer film structure, Stedim 40
- Pre-installed stirrer
- Ports for harvesting, sampling, probes, addition, aeration.
- Disposable sensors pre-installed
- Tubing and connectors: Application based / customized solutions
- Increase of safety via second exhaust filter connection
- Height/Diameter ratio (TV): 2:1
- Height/Diameter ratio (WV): 1.3:1

### 3. Impeller design Flexibility of mixing



- Type of impeller:
  - 2 x 3-blade segment impeller
  - 2 x 6-blade disk impeller
  - Combination
- Design is a copy of the classical impellers
- Impeller/Bag diameter ratio: 0.38

### 3. Impeller design Summary Mixing Studies 200L

liquid volume [l]	agitation [rpm]	mixing time [s]
2 x 3-blade segment impeller		
130	50	26
130	100	18
130	150	12
200	50	24
200	100	16
200	150	11

→ very efficient mixing!!!



### 3. Impeller design Summary Mixing Studies 200L

liquid volume [l]	agitation [rpm]	mixing time [s]
<b>2 x 6-blade disc impeller</b>		
130	50	20
130	100	15
130	150	13
200	50	25
200	100	18
200	150	16

→ very efficient mixing!!!


### 3. Impeller design 200L Tip speed & power input

agitation [rpm]	tip speed [m/s]
50	0.59
100	1.18
150	1.77

agitation [rpm]	power input per volume [W/m <sup>3</sup> ]	type of impeller
50	2.2	
100	17.5	
150	71.0	
50	4.0	
100	32.2	
150	132	

## 4. Aeration/gassing design

### Flexibility in gassing



#### 4. Aeration/gassing design Sparger design




μ-sparger (PE; 18-45 μm)




ring sparger (PC; 0.8 mm)


4. Aeration/gassing design :  
 Theoretical calculation  $k_La$  vs cell density (CHO cultivation example)



4. Aeration/gassing design : Summary kLa studies based on different impellers, stirrer speeds & flow rates




## Case study 1: CHO Cultivation



# Case Study 1

## CHO clone DG44 ST1-6 in Serum free chemically defined media


Viable cell density & Viability




⇒ No differences in cell density and viability between BIOSTAT® CultiBag STR and classical reusable bioreactor (typical 6-7 x 10<sup>6</sup> cells/mL for this clone & media)

# Case Study 2


## "Disposable Factory" Upstream





## Case study 3: CHO Cultivation (Customer X)



# Case Study 3 (Customer X) CHO Cultivation



- ◆ STR-4
- STR-5
- ▲ STR-6
- ✧ 2KL SS-4
- ✧ 2KL SS-5
- 2KL SS-6



⇒ Comparable Results

## Summary: „Classical“ single-use bioreactors



## Benefits for upscaling when utilizing single-use bioreactors with a classical design



- Scale-up parameters identical to classical process
  - tip speed, power input, agitation,  $k_L a$ ,...

- Easy optimization during the complete process
  - ⇒ User can implement his experience from classical processes into single-use applications



- Easy transfer to 5-10,000L stainless steel reactor due to comparability in design
  - ⇒ Reduction of time, need for optimization & cost during scale-up

## Future of modern pharmaceutical manufacturing for upstream processes

- Completely single-use stirred bioreactor
- Advanced control and data-storage
- Classical design on all critical parameters
  
- Easy transfer from reusable to single-use process
- Easy transfer from single-use to reusable process
- Scalability -> from R&D to production scale





Thank you for your attention.

Any Questions?

Dr. Andreas Kocourek, Sartorius Stedim Biotech, Application Specialist, Fermentation  
Australia ([andreas.kocourek@sartorius-stedim.com](mailto:andreas.kocourek@sartorius-stedim.com))

## 2. Bag Design Stedim 40

### 1 Polyethylene terephthalate (PET)

acts as a light, strong and clear protective outer layer. PET provides robustness..

### 2 Polyamide (PA)

increases durability and strengthens the bag..

### 3 Ethyl Vinyl Alcohol (EVOH)

acts as the main gas barrier.

### 4 Ultra-Low Density Polyethylene (ULDPE)

acts as the fluid contact layer. The S40 Polyethylene material is in compliance with usual pharmacopoeias and provides a clean, inert and highly chemical resistant contact layer.

